

Experiment # (1)

Measurements and Errors



Physical quantities:

- 1) Measured quantities (length, mass, and time).
- 2) Calculated quantities (volume, velocity, acceleration,)

Types of errors:

- 1) Personal error.
- 2) Systematic error (instrument error)
- 3) Random error.

*The error of measurement (Δx):

- 1) The smallest division for digital instruments.

Example: Balance: error = least count = 0.01g

- 2) Half smallest division for non digital instruments.

Example: Ruler, or meter stick: least count = 0.1 cm \rightarrow error = $\frac{0.1}{2}$ cm = 0.05 cm



1) if $z=x\pm y$ then: $\Delta z = \Delta x + \Delta y$, where x and y are measured quantities, z is calculated quantity(The errors are cumulative.)

2) if $z=xy$ or $z=x/y$, then:

a. take normal logarithm to both sides: $\ln(z)=\ln(x^*/y)=\ln(x) \pm \ln(y)$

b. derive explicitly : $\frac{\Delta z}{z} = \frac{\Delta x}{x} + \frac{\Delta y}{y}$

3) $z=a x^n$, where a , and n are constants, or unmeasured quantities.

$$\frac{\Delta z}{z} = n \frac{\Delta x}{x}$$

Any experimental value measured or calculated written as:
 $x\pm\Delta x$, where x is the value and Δx is the error in this value.


Example : Find the relative error in the volume of the cylinder?

$$V = \frac{\pi D^2 h}{4} \rightarrow \frac{\Delta V}{V} = \frac{2 \Delta D}{D} + \frac{\Delta h}{h}$$



*Graphics:

$y = a x + b$ where a : the slope, and b : y-intercept.



***Conclusions:**

I) Compare between real values and experimental values, why the similarities/differences?

II) Are you investigate the objectives of the experiment, how?

III) The errors that being in the experiment.

IV) The relation between x and y is linear and the constant of proportionality is




1) Vernier caliper:

For measure diameter of the cylinder (D).

$$\text{Least count} = \frac{1 \text{ mm}}{20} = 0.05 \text{ mm}$$

$$\text{Error} = \frac{0.05}{2} \text{ mm} = 0.025 \text{ mm}$$



The higher valued of the vernier scale mark lines up with the main scale is 7.

The reading of the Vernier caliper = 1.1cm + 7 * 0.1 mm = 11mm + 0.7 mm = 11.7 mm.




II) Micrometer: (most accurate instrument to measure length)

For measure height of the cylinder (h).

$$\text{Least count} = \frac{0.5 \text{ mm}}{50} = 0.01 \text{ mm}$$

$$\text{Error} = \frac{0.01}{2} \text{ mm} = 0.005 \text{ mm}$$




The reading of the Micrometer = $5.5 \text{ mm} + 27 \times 0.01 \text{ mm} = 5.77 \text{ mm}$.




The density of the cylinder:

Density = $\frac{\text{mass}}{\text{Volume}}$ $\rightarrow \rho = \frac{m}{\frac{\pi D^2 h}{4}}$



#	h (mm) ±0.005	D (mm) ±0.025	V(cm ³)	m (g) ±0.01
1				
2				
3				
4				



THANK YOU

